



GETTING READY For Your CGM Appointment

Congratulations on beginning your journey with the Eversense 365 CGM System.
Let's get you started!

BEFORE Your Appointment

Download the Eversense 365 Mobile App & create an account

- ☐ From your compatible mobile device, download the Eversense 365 App from the App Store® or Google Play™

• Download this app → 

- ☐ Follow the prompts to create your Eversense account
 - On the Welcome screen, choose "I do not have a Smart Transmitter"
 - On the Unit of Measurement screen, leave it set at mg/dL
 - Write down your account username (email) and password:
email: _____
password: _____
- ☐ After creating your account, you will receive an email from Eversense. Click on the link in email within 72 hours to activate your account.

Questions? We're here for you!
Eversense Customer Support:
1-844-736-7348 (1-844-SENSE4U)

DAY OF Your Appointment

- ☐ Wear a loose, short-sleeved or sleeveless shirt
- ☐ What to bring:
 - Your mobile device, fully charged, with the Eversense 365 App installed.
 - Log into the app with your Eversense account username (email) and password.
- ☐ You can return to your normal activities. However, to allow proper healing, it's not recommended that you soak in a tub, swim, or do strenuous activities for five days.
- ☐ Make sure to have a blood glucose meter, valid test strips, and lancets at home so that you can begin calibrating the system after warm up is complete.

Check out our website for videos on setting up your CGM System:
www.eversensecg.com/resources-365



Sign up for your live Eversense online session today!
Our clinical team will prepare you for sensor placement and get you off to a great start using your Eversense 365 system. Don't miss out!

◀ **Scan here to view upcoming dates and secure your spots before your appointment.**



© 2026 Senseonics, Incorporated.
All rights reserved.

Manufactured by
Senseonics, Incorporated.
20451 Seneca Meadows Parkway
Germantown, MD 20876-7005 USA
eversensecg.com

The Eversense® 365 Continuous Glucose Monitoring (CGM) System is indicated for continually measuring glucose levels for up to one year in people (18 years and older) with diabetes. The system is indicated for use to replace fingerstick blood glucose (BG) measurements for diabetes treatment decisions. Fingerstick BG measurements are required for calibration one time a week after day 13, and when symptoms do not match CGM information or when taking medications of the tetracycline class. The sensor insertion and removal procedures are performed by a health care provider. The Eversense 365 CGM System is a prescription device; patients should talk to their health care provider to learn more.

For safety information, see bit.ly/eversensesafety

Eversense, Eversense 365 Continuous Glucose Monitoring, and the Eversense logo are trademarks of Senseonics, Incorporated. Other brands and their products are trademarks or registered trademarks of their respective holders and are used solely for informative purposes. No relationship or endorsement should be inferred or implied. MKT-001725 Rev 1