

Welcome to
**ONE YEAR.
ONE CGM.**



Congratulations on Getting Eversense 365 CGM!

TO SUPPORT PROPER HEALING, FOLLOW THESE STEPS OVER THE NEXT FIVE DAYS:



Do not swim or soak in a tub for at least five days. Showering is fine but avoid the direct stream of water on the site.



The Tegaderm™ should stay on for at least five days. If it comes off before then, replace it. Take extra care not to pull off the Steri-Strips™ when removing the Tegaderm™ pad.



Until the site heals, **avoid activities that may pull at the incision site or cause sweating** around the insertion area.



Leave Steri-Strips™ on until they fall off. Trim the edges of the Steri-Strips™ if they start to curl; do not remove them when doing so.

NOTIFY YOUR HEALTHCARE PROVIDER IF:

Your Steri-Strip™ comes off before your incision is fully closed or if you develop a fever or experience pain, redness, swelling, warmth or drainage at the incision site.

Step 1: Warm-Up

- After the sensor is linked with the smart transmitter, the 24-hour warm-up phase begins.
- You do not need to wear your smart transmitter during warm-up.
- Charge your smart transmitter during this time so that it is fully charged and ready to use.

Step 2: Initialization Phase

- Once warm-up is complete, place the smart transmitter over the sensor with an adhesive patch using the Placement Guide.
- Complete 4 calibrations, 2 to 12 hours apart. Refer to the Quick Reference Guide on how to enter a fingerstick value.
- After the 2nd successful calibration, you will begin to see your glucose data.
- Enter your CGM settings provided by your health care provider into the mobile app.

Step 3: Daily Wear

- After you enter the 4th calibration, Daily Calibration phase begins. Calibrate once daily for 13 days.
- At Day 14, you just need to calibrate one time per week for one year of accurate CGM results!
- You'll need to charge the transmitter and replace the adhesive patch every day. Try to do it the same time every day, like while taking your morning shower.

Who to Contact with Questions

- For concerns regarding your insertion site, contact the healthcare provider who did your procedure.
- For diabetes management questions, contact the healthcare provider who prescribed Eversense.
- For technical questions or help getting started with Eversense, contact **1-844-SENSE4U (736-7348)**.



Find helpful videos are available at: eversensecgm.com/resources